

| | | |
|---|----------------------------------|--|
|  | Product Information Sheet | Page 1 of 3 |
| | PIS-060 EB-TRK-Preps | Version: 002 Effective Date: 07/06/20 Author: Beth Lowe CO#: 062520-1 |

Print only page 2-3 for Customers

Concentration: Not Applicable

Storage and Handling:

Store at room temperature upon arrival. During shipment, crystals or precipitation may form on the TRK Buffer. Dissolve by warming buffer to 37°C.

Ordering Information:

| Item Number | Preparations | Number of Tubes and Volumes |
|-------------|--------------|---|
| EB-TRK-50 | 50 | EB-TRK DEPC Water: 1x10mL, EB-TRK RNA Wash Buffer I: 1x50mL, EB-TRK RNA Wash Buffer II: 1x12mL, EB-TRK TRK Lysis Buffer: 1x40mL, EB-TRK 2mL Collection Tubes: 100, EB-TRK RNA Mini Columns:50 |
| EB-TRK-200 | 200 | EB-TRK DEPC Water: 1x40mL, EB-TRK RNA Wash Buffer I: 1x200mL, EB-TRK RNA Wash Buffer II: 1x50mL, EB-TRK TRK Lysis Buffer: 1x150mL, EB-TRK 2mL Collection Tubes: 400, EB-TRK RNA Mini Columns: 200 |

Product Description:

EB Pure Total RNA Kit is a convenient system for the fast and reliable purification of RNA from cells or soft tissue. The key to this system is the High-Bind matrix that specifically, but reversibly, binds RNA under optimized conditions allowing proteins and other contaminants to be removed. Then RNA is easily eluted with deionized water or a low salt buffer. Each column can bind up to 100µg of RNA. This kit is used to recover RNA ready for applications such as RT-PCR, Northern blotting, poly A+ RNA (mRNA) purification, nuclease protection, and in vitro translation.

Additional Materials and Equipment to be Supplied by User:

- Microcentrifuge capable of at least 14,000 x g
- RNase-free pipettes and 1.5 mL microcentrifuge tubes
- 100% ethanol
- 70% ethanol in sterile DEPC-treated water
- Homogenization equipment
 - Homogenizer columns
 - Needle and Syringe
 - Mortar and pestle
 - Glass beads
 - Rotor-stator Homogenizer
- Optional: 14.3M 2-mercaptoethanol (β-mercaptoethanol)

Preparation Procedure:

Prepare RNA Wash Buffer II according to the table below and store at room temperature:

- For EB-TRK-50 Kits add 48mL of 100% Ethanol to RNA Wash Buffer II Bottle and indicate addition on the bottle.
- For EB-TRK-200 Kits add 200mL of 100% Ethanol to each RNA Wash Buffer II Bottle and indicate addition on the bottle.
- Optional: Add 20 µL 2-mercaptoethanol (β-mercaptoethanol) per 1 mL TRK Lysis Buffer. This mixture can be stored for 2 weeks at room temperature.

Procedure Options:

Centrifugation Protocol – Listed below and in Lab Manual available online at empiricalbioscience.com

Vacuum Protocol – Listed in the Lab Manual available online at empiricalbioscience.com

*This product is intended for Research Use Only. This product is manufactured under ISO13485:2016 Quality System Requirements and is available for use as a Raw Material for use in IVD applications. Please contact Empirical Bioscience for further details.

For MSDS and Certificate of Analysis please visit www.empiricalbioscience.com

Centrifugation Protocol

- 1. **For Tissue:** Determine the proper amount of starting material.
 - Homogenize and disrupt the tissue.
 - Add 350 μL of TRK Lysis Buffer for less than 15 mg of tissue, or 700 μL for 20-30 mg of tissue.
 - Centrifuge and maximum speed ($\geq 12,000 \times g$) for 5 minutes.
 - Transfer only the clarified supernatant to a fresh microcentrifuge tube. Avoid any fatty upper layer that may have formed.**For Cultured Cells:**
 - Determine the proper amount of starting material.
 - Harvest the cells and resuspend the pellet in 350 μL of TRK Lysis Buffer for less than 5×10^6 cells, or 700 μL for $5-10 \times 10^6$ cells.
 - Homogenize the cells using a syringe and needle or homogenizer mini-columns.**Optional:** Add 20 μL β -ME per 1mL TRK Lysis Buffer before using.
- 2. Add 1 volume 70% ethanol. Vortex to mix thoroughly. Do not centrifuge. If the sample has lost any volume during homogenization, adjust the volume of ethanol to match.
- 3. Insert a High-Bind RNA Mini Column into a 2mL Collection Tube (provided).
- 4. Transfer 700 μL sample to the Mini Column. Centrifuge at 10,000 $\times g$ for 1 minute. Discard the filtrate and reuse the Collection Tube.
- 5. Repeat Step 4 until all of the sample has been transferred to the column.
- 6. **Optional:** Perform DNase digestion according to your preferred protocol.
- 7. Add 500 μL RNA Wash Buffer I to the column. Centrifuge 10,000 $\times g$ for 30 seconds. Discard the filtrate and reuse the Collection Tube.
- 8. Add 500 μL RNA Wash Buffer II diluted with 100% ethanol (See Bottle) to the column. Centrifuge 10,000 $\times g$ for 1 minute. Discard the filtrate and reuse the Collection Tube.
- 9. Repeat Step 8.
- 10. Centrifuge at maximum speed for 2 minutes to dry the column.
- 11. Transfer the High-Bind Column to a clean 1.5mL microcentrifuge tube (not provided).
- 12. Add 40-70 μL DEPC Water. Centrifuge at maximum speed for 2 minutes.
- 13. Store eluted RNA at -70°C .

***This product is intended for Research Use Only. This product is manufactured under ISO13485:2016 Quality System Requirements and is available for use as a Raw Material for use in IVD applications. Please contact Empirical Bioscience for further details.**

For MSDS and Certificate of Analysis please visit www.empiricalbioscience.com