

	Product Information Sheet	Page 1 of 2
	PIS-056 MgCl₂-25-Volume	Version: 001 Effective Date: 07/06/20 Author: Beth Lowe CO#: 062520-1

Print Page 2 Only

Product Information Sheet

Product Name: 25mM Magnesium Chloride*

Concentration: Magnesium Chloride: 25mM

Storage and Handling:

Upon arrival store at -20°C for provided expiration date. This product can be stored at Room Temperature for 120 days, 4°C for 240 days.

Ordering Information:

Item Number	Number of Tubes and Volume
MgCl2-25-25mL	25mM Magnesium Chloride: 25x1mL

Product Description:

The Magnesium Chloride is supplied as a 25 mM solution. This product is free of any DNase or RNase activity and is also free any contaminating E.coli or Human genomic DNA.

*This product is intended for Research Use Only. This product is manufactured under ISO13485:2016 Quality System Requirements and is available for use as a Raw Material for use in IVD applications. Please contact Empirical Bioscience for further details.
For MSDS and Certificate of Analysis please visit www.empiricalbioscience.com