

Agarose Gel Extraction Kit

Enzymes and Reagents

PCR | qPCR | RT | Kits | Oligonucleotides



Empirical Bioscience Agarose Gel Extraction Kit is designed for high-yield recovery of DNA from agarose gels with simultaneous removal of primer-dimers, primers, nucleotides, proteins, salt, agarose, ethidium bromide, and other impurities.

The preparations are based on silica-membrane technology for binding DNA in high-salt and elution in low-salt buffer. The kit provides a simple and efficient way to purify DNA in a size range between 100 bp and 10 kb. The extraction kit requires no organic extraction or precipitation and guarantees high and stable recovery rates.

Item Number	Preparations
SS-PPK-10	10
PPK-50	50
PPK-250	250

Features

Wide Range

- Purify DNA fragments from 100bp to 10kb.

Easy-to-Use

- Spin Columns & Tubes.
- No organic extractions or precipitations necessary.

Robust Performance

- High yield recovery of DNA.

Kit Components:

- Extraction Buffer
- Activation Buffer
- Washing Buffer
- Elution Buffer
- Spin Columns & Collection Tubes

Related Products:

- Taq DNA Polymerase
- Taq 2X Master Mix
- Taq 2X MeanGreen Master Mix
- FlashTaq HotStart DNA Polymerase
- FlashTaq HotStart 2X Master Mix
- FlashTaq HotStart 2X MeanGreen Master Mix
- 5X MeanGreen Buffer
- EasyLoad Low Range DNA Ladder
- EasyLoad Mid Range DNA ladder
- EasyLoad High Range DNA Ladder
- UltraPure dNTP Solution Mix, 10mM each
- UltraPure dNTP Solution Set, 100mM each
- 6X MeanGreen Loading Dye



Observe the Difference

TO ORDER, REQUEST SAMPLES OR FOR MORE INFORMATION, CONTACT:

Empirical Bioscience, Inc. | 2007 Eastcastle Drive SE

Grand Rapids, MI 49508 | Phone: 877.479.9949 Fax: 616.825.6149

Email: sales@empiricalbioscience.com

empiricalbioscience.com

Empirical Bioscience Product Line

Taq DNA Polymerases	
ITEM NO.	DESCRIPTION
TP-NB-100	Taq DNA Polymerase without Buffer, 100 Units
TP-NB-500	
TP-NB-1000	
TP-NM-100	Taq DNA Polymerase with 10X Reaction (no MgCl ₂) Buffer
TP-NM-500	
TP-NM-1000	
TP-MG-100	Taq DNA Polymerase with 10X Reaction (MgCl ₂) Buffer
TP-MG-500	
TP-MG-1000	

FlashTaq™ DNA Polymerases	
ITEM NO.	DESCRIPTION
FT-NB-100	FlashTaq™ DNA Polymerase without Buffer
FT-NB-500	
FT-NB-1000	
FT-NM-100	FlashTaq™ DNA Polymerase with 10X Reaction (no MgCl ₂) Buffer
FT-NM-500	
FT-NM-1000	
FT-MG-100	FlashTaq™ DNA Polymerase with 10X Reaction (MgCl ₂) Buffer
FT-MG-500	
FT-MG-1000	

Integrity DNA Polymerases™	
ITEM NO.	DESCRIPTION
HF-RB-MG-100	Integrity High-Fidelity DNA Polymerase™ with 10X Reaction Buffer
HF-RB-MG-200	
HF-RB-MG-500	
HF-RB-MG-1000	

PFU-50 High Fidelity Polymerase	
ITEM NO.	DESCRIPTION
PFU-50-100	PFU-50 High Fidelity Polymerase (50X Taq) with Reaction Buffer
PFU-50-200	
PFU-50-500	
PFU-50-1000	
PFU-50-5000	

dNTP Solutions	
ITEM NO.	DESCRIPTION
dNTP-1mL	UltraPure Deoxynucleotide Solution Mix, 10mM
dNTP-100-S0.2mL	UltraPure Deoxynucleotide Solution Set, 100mM
dNTP-100-S1mL	
dNTP-100-S5mL	

Taq Master Mixes	
ITEM NO.	DESCRIPTION
TP-MM-100	Taq 2X Master Mix
TP-MM-250	
TP-MM-500	
TP-MM-1000	
TP-MM-2500	
TP-MMWD-100	Taq 2X MeanGreen™ Master Mix
TP-MMWD-250	
TP-MMWD-500	
TP-MMWD-1000	
TP-MMWD-2500	5X Multiplex Master Mix
TP-MP-100	
TP-MP-500	
TP-MP-1000	
TP-MP-GCE-100	Multiplex Master Mix with GC Enhancer
TP-MP-GCE-500	

FlashTaq™ Master Mixes	
ITEM NO.	DESCRIPTION
FT-MM-100	FlashTaq™ 2X Master Mix
FT-MM-250	
FT-MM-500	
FT-MM-1000	
FT-MMWD-100	FlashTaq™ 2X MeanGreen™ Master Mix
FT-MMWD-250	
FT-MMWD-500	
FT-MMWD-1000	

EasyLoadDNA Ladders™	
ITEM NO.	DESCRIPTION
LA-EL-LR-100	EasyLoad™ Low-Range DNA Ladder, 100 Lanes
LA-EL-MR-100	EasyLoad™ Mid-Range DNA Ladder, 100 Lanes
LA-EL-HR-100	EasyLoad™ High-Range DNA Ladder, 100 Lanes

AlloyTaq Blends™	
ITEM NO.	DESCRIPTION
AT-HF-RB-100	AlloyTaq™ High-Fidelity Blend with 10X Reaction Buffer
AT-HF-RB-200	
AT-HF-RB-500	
AT-HF-RB-1000	

GC Enhancers	
ITEM NO.	DESCRIPTION
GCE-10-100	10X GC Enhancer
GCE-10-500	

MeanGreen™ Buffer	
ITEM NO.	DESCRIPTION
RB-WD-20mL	5X MeanGreen™ Buffer (20 x 1mL)

MeanGreen™ Loading Dye	
ITEM NO.	DESCRIPTION
LD-MG-4mL	6X MeanGreen™ Loading Dye

rEVALution™ qPCR Mixes	
ITEM NO.	DESCRIPTION
RQ-100	rEVALution 2X™ qPCR Master Mix with separate ROX®
RQ-500	
RQ-1000	
RQ-2500	

2x qPCR Probe Master Mix™	
ITEM NO.	DESCRIPTION
QP-MM-100	2x qPCR Probe Master Mix with separate ROX®
QP-MM-500	
QP-MM-1000	
QP-MM-2500	

RTScript™	
ITEM NO.	DESCRIPTION
RT-CSK-100	Reverse Transcriptase cDNA Synthesis Kit

Kits	
ITEM NO.	DESCRIPTION
PMK-50	Plasmid Mini-Prep Kit
PMK-250	
PPK-50	PCR Purification Kit
PPK-250	
AGEK-50	Agarose Gel Extraction Kit
AGEK-250	

Oligonucleotides	
ITEM NO.	DESCRIPTION
RTG-PR	Ready to Go™ Primers