

# EasyLoad™ DNA Ladders




Enzymes and Reagents

PCR | qPCR | RT | Kits | Oligonucleotides

All Empirical EasyLoad™ DNA Ladders contain MeanGreen™ loading dye.

	Low Range	Mid Range	High Range
<b>Range</b>	50bp-1kb	100bp-3kb	500bp-10kb
<b>Band Sizes</b>	50, 75, 100, 150, 200, 300 400, 500, 600, 800, 1000bp	0.1, 0.15, 0.2, 0.3, 0.4, 0.5, 0.6, 0.8, 1, 1.5, 2, 3kb	0.5, 0.6, 0.8, 1, 1.5, 2, 3, 4 5, 6, 8, 10kb
<b>Conc.</b>	120ng/μl	130ng/μl	120ng/μl
<b>Dye Sizes</b>	Blue ~ 4000bp Yellow <25bp	Blue ~ 4000bp Yellow <25bp	Blue ~ 4000bp Yellow <25bp



<b>Lanes/Volume</b>	100 Lanes/500μl	100 Lanes/500μL	100 Lanes/500μL
<b>Item #</b>	<b>LA-EL-LR-100</b>	<b>LA-EL-MR-100</b>	<b>LA-EL-HR-100</b>



**TO ORDER, REQUEST SAMPLES OR FOR MORE INFORMATION, CONTACT:**  
 Empirical Bioscience, Inc. | 2007 Eastcastle Drive SE  
 Grand Rapids, MI 49508 | Phone/Fax: 877.479.9949  
 Email: [sales@empiricalbioscience.com](mailto:sales@empiricalbioscience.com)  
[empiricalbioscience.com](http://empiricalbioscience.com)