



## Product Information Sheet

**Product Name:** Deoxynucleotide Solution Mix, 10mM each\*

**Item No.:** dNTP-Volume, SS-dNTP-Volume

**Concentration:** 10mM each

**Storage and Handling:**

Store at -20°C upon arrival.

**Purity:**

≥ 99% (HPLC)

**Product Description:**

Empirical's Deoxynucleotide Solution Mix is an equimolar solution of ultrapure dATP, dCTP, dTTP, and dGTP supplied as clear aqueous solution (pH 8.5).

**Applications:**

For standard PCR applications a final concentration of 200µM each dNTP is recommended.

**\* This product is for “Research Use Only. Not for use in diagnostic procedures”.**

For MSDS and Certificate of Analysis please visit [www.empiricalbioscience.com](http://www.empiricalbioscience.com)