

25mM Magnesium Chloride

Enzymes and Reagents

PCR | qPCR | RT | Kits | Oligonucleotides

The Magnesium Chloride is supplied as a 25 mM solution. This product is free of any DNase or RNase activity and is also free any contaminating E.coli or Human Genomic DNA.

Storage and Handling:

Recommended to store at 20°C for long-term storage. This product can be stored at 2°C - 8°C for short-term storage.

Item Number	Total Volume Received	Quantity Received
SS-MgCl ₂ -25-0.25mL	250 µl	(1x250µL)
MgCl ₂ -25-25mL	25ml	(25x1 mL)